INSTRUCTION MANUAL
ELECTRONIC BATHROOM SCALE

B18- / BS-108

Distributed by; V.MAX North America
Manufactured in; PRC Mainland

YOUR NEW ELECTRONIC PERSONAL SCALES HAVE BEEN DESIGNED TO ACCURATELY INDICATE YOUR WEIGHT GAIN OR LOSS OVER A PERIOD OF TIME AND SHOULD GIVE YOU MANY YEARS OF SERVICE UNDER NORMAL USE

Specifications

<ul style="list-style-type: none">• Anti-Slippery lateral strips• Equipped with high precision "strain-Gauge" sensors• Capacity 180kg/396lb/28.8st• Division 100g/0.2lb• Backlight	<ul style="list-style-type: none">• LCD Display• Automatic Zero Resetting• Automatic Switch Off• Over Load Indication• 2 batteries CR2032
--	---

Power Supply

Please insert batteries or remove the insulating strip from the battery contacts before using. Replace the batteries when the "Low Battery Indicator" shows : 

Operating Instructions

1. Place the scale on a hard, flat surface, this will ensure the greatest accuracy and repeatability.
2. Stepping onto the scale, then the number on LCD flash 3 times, that is the weighing result of your body.

Automatic Switch Off

If no weighing is carried out while the zero is being displayed, the scale will automatically switch off after approximately 10 seconds.

Overload Indication

When the scales are overloaded (over 180kg/396lb/28st), the following will appear on the screen.

Battery Change Indication

When the batteries need replacing, the following will appear on the screen:

Please open the back board of the scales, then take batteries carefully out.

Care and Attention

- Place your scales on a hard, flat surface, this will ensure the greatest accuracy and repeatability. Uneven floors may affect the readings.
 - Always weigh yourself on the same scale placed on the same floor surface.
 - Protect the scale against shocks, damp, dust, chemicals, drastic changes in temperature, electromagnetic field and source of heat that are too near.
 - Never step onto the outer edge of the scale with one foot, the scale could flip.
 - Keep packaging material away from children!
 - Do Not dispose of the old batteries in fire. The batteries may explode or leak. Recycle or manage as hazardous waste.
 - Do not mix Old Batteries and New Battery. Do not mix Alkaline, Standard (Carbon-Zinc) or rechargeable (Nickel-Cadmium) Battery.
 - If the scale would not be used in a long period, please take the battery out.
 - **Cleaning:** You can clean the scale with a damp cloth and a little washing liquid, if required. Never immerse the scale in water or rinse it under running water.
-

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